



### References

- Archer Daniels Midland
- Cardinal Stritch College
- Cambridge Major Laboratories
- Cedarburg Pharmaceuticals
- Columbus Chemical Industries, Inc.
- Eastman Chemical Company
- Eli Lilly & Company
- Epicentre
- Gala Design, Inc.
- Harlan Bioproducts for Science
- Hydrite Chemical Company
- Kraft
- Laboratory Corporation of America
- Lucigen Corporation
- Madison Area Technical College
- Medical College of Wisconsin
- Milwaukee Metropolitan Sewerage District
- Monsanto
- Northern Illinois Police Crime Laboratory
- Nutrasweet Company
- Pierce Biotechnology
- S.C. Johnson & Sons, Inc.
- Safety Kleen Corporation
- Scientific Protein Labs
- Stepan Company
- Wisconsin Veterinary Diagnostic Laboratory

### Consulting

Ironically, companies and universities don't want to experiment with laboratories. That's why so many turn to LabTech Midwest for our expert consulting services.

We're up to any challenge. We can help you define your project, design it, and get what you want most in a laboratory — a guaranteed result.

### Close Enough for Comfort

From our headquarters just west of Milwaukee, in New Berlin, WI, we serve customers throughout the Midwest, especially across Southeastern Wisconsin and Northeastern Illinois.

In other words, we're close enough to be there when you need us. Over the years, our combination of world-class expertise and local service has been why so many companies have trusted Labtech Midwest with testing, certifying and servicing thousands of laboratory containment devices.

To find out more, phone, fax or e-mail us. Or visit our web site.



2100 S. Calhoun, New Berlin, WI 53151  
 Phone: 262-754-9800 Fax: 262-754-9802  
 e-mail: mbishop@labtech-midwest.com www.labtech-midwest.com

**TESTING & CERTIFICATION**  
**LABORATORY DESIGN • SERVICE & REPAIR • CONSULTING**

**The DANGER is invisible,  
odorless, silent and intangible.**



**But it  
can't hide  
from LabTech  
Midwest.**





## LabTech Midwest, Inc. The Safe Lab People.

Nobody can give you peace of mind about lab safety like LabTech Midwest, Inc. For more than ten years, we've been testing and servicing fume hoods, biological safety cabinets and other laboratory containment devices.

Our testing is in full compliance with the latest, most rigorous standards in the industry. (Astonishingly, some companies are still testing under guidelines five years old, or more.)

## Trust the People the Wisconsin Biotechnology and Medical Device Association Trusts.

It's one thing for us to say how good we are. But it's quite another for the Wisconsin Biotechnology and Medical Device Association to agree with us. They recently chose LabTech Midwest to provide testing and certification services for their members.

We take this as a testament to our quality, thoroughness and commitment to excellence.

## A Full Range of Services.

In more than a decade, our company has grown and our services have greatly expanded. We offer Testing & Certification, of course. But also Service and Repair, Laboratory Design, and Consulting.

## Testing & Certification

We work to the highest standards. We are trained and equipped to meet and exceed the demands of OSHA and the NEW Laboratory ventilation Standard (ANSI Z9.5). Our comprehensive Testing and Certification services for all laboratory containment devices include:

- Fume Hoods — ASHRAE 110, SEFA 1, ANSI Z9.5 and LTMW 1104 SOP.
- Biological Safety Cabinets — NSF Standard 49
- Laminar Flow Device Certification — Federal Standard 209B and IEST-RP-CC0002.2
- Clean Rooms — ISO 14644, and Federal Standard 209E — NEEB trained.
- Other Testing — We will work with you to develop custom testing methodology for any application.

## Service & Repair

We are committed to delivering complete care for today's modern laboratory. No project is too large — or too small. We offer complete service and repair on fume hoods, biological safety cabinets, laminar flow devices, glove boxes, incubators and more.

Many of our best customer relationships began with providing great service.

## Laboratory Design & Testing

We can be your single-source solution to lab planning, design, products and installation ... all the way to final move in. We offer a complete turnkey service.

Plus, we can save you money. In addition to new products, we buy and sell gently used fume hoods and casework. This can often give you the highest quality at the best price.

